Manuale di Istruzione Uso e Manutenzione
User and Maintenance manual
Manuel d'Utilisation et d'Entretien
Gebrauchs- und Wartungsanleitung
Manual de Instrucciones de Uso y Mantenimiento

Sterilintraplant Plus

carlo de giorgi s.r.l.

Dati modello**Model data / Donnés du modèle / Modelldaten / Datos del modelo**

Modello Model Modèle Modell Modelo	Alimentazione Power supply Alimentation Betriebs-spannung Alimentación	Anno di costruzione Year of manufacture Année de fabrication Baujahr Año de fabricación
557/00 558/00	<input type="checkbox"/> 230 V 50 Hz <input type="checkbox"/> 120 V 60 Hz	<input type="checkbox"/> 2004 <input type="checkbox"/> 2005 <input type="checkbox"/> 2006 <input type="checkbox"/> 2007 <input type="checkbox"/> 2008

**carlo de giorgi** s.r.l.

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Tel. +39.02.356.15.43 – Fax +39.02.356.18.08

CONTENTS

1 - STANDARDS AND GENERAL PRECAUTIONS	3
1.1 - INTRODUCTION	3
1.2 - DESCRIPTION OF THE APPLIANCE	3
1.3 - REFERENCE STANDARDS	3
1.3.1 - Obligatory compliance	3
1.3.2 - Voluntary compliance	4
1.4 - RESPONSIBILITIES OF THE CUSTOMER.....	4
1.5 - REQUESTS FOR SERVICE AND SPARE PARTS.....	4
2- TECHNICAL FEATURES.....	4
3 - COMPONENTS.....	4
4 - OPERATION AND USE	5
4.1 - OPERATION	5
4.2 - OPERATOR	8
4.3 – INTENDED USE	8
4.4 - SAFETY WARNINGS.....	8
5 - OPERATOR INSTRUCTIONS	9
5.1 - CONTROLS.....	9
5.2 - INSTRUCTIONS FOR USE	11
5.2.1 - Setting into operation.....	11
5.2.2 - Normal deactivation and emergency deactivation procedures.....	11
5.2.3 – Editing a program	11
6 - MAINTENANCE	12
6.1 - CLEANING	12
6.2 - STERILIZATION.....	13
7 - DIAGNOSTICS.....	13
8 - ACCESSORIES.....	14
8.1 - MOTOR	14
8.2 – TUBE FOR IRRIGATION	14
8.3 - UNION	14

1 - STANDARDS AND GENERAL PRECAUTIONS

1.1 - INTRODUCTION

This manual is the property of CARLO DE GIORGI S.r.l. The contents of this manual cannot be reproduced or disclosed to other parties. All rights reserved.

1.2 - DESCRIPTION OF THE APPLIANCE

The appliance is a micromotor for use in the dental sector for surgery and implantology, only when equipped with an implantology handpiece.

This device must be used exclusively by suitably trained persons who are authorised to perform the practices in question and who have accumulated sufficient experience in the use of appliances of this type.

Before using the micromotor it is mandatory to read the contents of this manual in their entirety, learning about the position of and function performed by all the controls and paying due attention to all the warnings and precautions that are prescribed. Failure to observe the above obligatory instruction or improper use of the micromotor may cause injury to the operator, to the patient, or to other persons present in the area.

IT IS STRICTLY PROHIBITED TO ATTEMPT TO OPEN OR REPAIR THE APPLIANCE; IN THE EVENT OF FAULTS OR ANOMALIES ALWAYS CONTACT CARLO DE GIORGI.

1.3 - REFERENCE STANDARDS

1.3.1 - Obligatory compliance

- EEC Directive n°.93/42 – Medical Devices (Italian Legislative Decree D.Lgs. n° 46/1997)
- EEC Directive n°.98/37 – Machinery Directive (Italian Presidential Decree D.P.R. n° 459/1996)
- EEC directive n°.73/23 – Low voltage Directive (Italian legislation: Law n° 791/1977, Legislative decrees D.Lgs n° 626/1996 and D.Lgs. n° 277/97)
- EEC directive n° 89/336 relative to Electromagnetic Compatibility (EMC) - (Italian legislation: Legislative decree D.Lgs n° 615/1996)
- DPR 27.4.1955 n°.547, "Norme per la prevenzione degli infortuni sul lavoro" (work accident prevention regulations)
- DPR 27.4.1956 n°.303, "Norme generali per l'igiene del lavoro" (general work hygiene regulations)
- Italian Legislative Decree 12 August 1991 n°.277 implementing EEC directives n°.80/605; n°.83/477; n°.86/188; n°.88/642 concerning the protection of workers against risks deriving from exposure to chemical, physical and biological agents in the workplace, as defined by art. 7 of Italian Law 30 July 1990, n°.212.

1.3.2 - Voluntary compliance

- EN 292 (1992) Machine safety – Basic concepts; general design principles – Part 1 – Basic terminology and methods (EN 292-1) /Part 2 – Specifications and technical principles (EN 292-2)

1.4 - RESPONSIBILITIES OF THE CUSTOMER

The user must install the appliance in a suitable place equipped with an electrical system that complies with established legislation.

The appliance must be installed in a dry environment that is illuminated in accordance with the requirements of established legislation.

NOTE: The expression "established legislation" refers to the laws in force in the user's country.

1.5 - REQUESTS FOR SERVICE AND SPARE PARTS

For all electrical maintenance requirements, contact CARLO DE GIORGI of Baranzate di Bollate, Via Tonale, 1.

In the case of faults, notify CARLO DE GIORGI technical personnel of the exact nature of the problem, specifying also the lot number shown on the label affixed to the rear of the appliance.

2- TECHNICAL FEATURES

		557/00 558/00
Voltage	V	230/120
Frequency	Hz	50/60
Max. power consumption	W	120
Nominal output voltage	V	30
Delayed acting fuse T	A	2x1.6 (230V)
Mass of the appliance	kg	6

Class II app.

type BF

The appliance features double electrical insulation (Class II).

3 - COMPONENTS

The package contains:

- a) power supply unit;
- b) power cord;
- c) micromotor (*);
- d) support for micromotor (*);
- e) footswitch;

- f) disposable irrigator;
- g) saline solution pipeline connection;
- h) stand for saline solution flask;
- i) flask of saline solution with plastic cradle;
- j) n° 9 wrench;
- k) user instructions;
- l) motor endcap (*);

(*) 1 on model 557/00, 2 on model 558/00.

The appliance can be moved by hand without requiring lifting devices.

4 - OPERATION AND USE

4.1 - OPERATION

The power supply unit must be placed on a surface next to the operator in such a way that:

- it is in a stable position without risk of falling;
- the micromotor connection lines are not subjected to excessive strain – risk of breakage of the wires or dragging of the unit;
- the functional status of the appliance can be checked easily.

Referring to Fig. 1, insert micromotor supports (1) in the lateral sockets; connect the micromotor to the socket (2-3) by inserting the plug and then twisting the locking ring without exerting undue force;

Fig. 1

With reference to Fig. 2, connect the footswitch to socket (4) on the rear panel by inserting the plug and twisting the locking ring; Insert the saline solution stand into upper socket (8) and screw it down, securing it by tightening the nut (using the supplied wrench) in the required position; now fit the saline solution bottle inside in the plastic cradle.


Fig. 2

With reference to Fig. 3 open the side door of peristaltic pump (1) and remove the tube (keeping it safe for future uses with irrigators other than the one supplied); insert the larger section tube (2) (included with the disposable irrigator supplied) into the peristaltic pump, paying attention to:

- the positioning of supports (3), since, as can be observed in the figure, the supports perform a restraining function on the exterior of the pump;
- the connection of the "IN" and "OUT" sides (4) respectively to the saline solution bottle and the surgical handpiece (not supplied).






Fig. 3

If the handpiece is of the dual irrigation type (see Fig. 4), connect the double outlet union (Fig. 5-1) to the irrigator supplied; in the event of single irrigation types, connect the single outlet union (Fig. 5-2). Fig. 5 shows the two types of unions supplied with the irrigator in the package.

**Fig. 4****Fig. 5****Fig. 5-1****Fig. 5-2**

Finally, plug the power plug into the mains socket outlet, connecting the appliance side (see Fig. 2) to socket (5) located on the rear of the unit; the same figure shows the main power switch (6) and the two fuses (7). At this point attach the required surgical handpiece (not included) to the micromotor.

Always switch OFF the appliance after use by setting the main power switch to "O".

4.2 - OPERATOR

After selecting the required program the operator starts the appliance by means of voluntary activation of the footswitch thereby assuming direct control of its operation.

NOTE: For production reasons the 32 preset programs contain standard parameters that are not of specific relevance for the implantologist; the programs must therefore be personalised in relation to individual requirements.

4.3 – INTENDED USE

The implantological micromotor is designed for use in dental and veterinary surgeries for surgical and implantological operations.

The appliance is not designed for any uses other than those specified by the manufacturer.

4.4 - SAFETY WARNINGS

- The appliance must be used exclusively by skilled personnel who are aware of the potential hazards;
- Do not use the appliance in the vicinity of inflammable materials or gases;
- Do not open or disassemble the appliance;
- Wear suitable work clothing and observe the safety regulations in force in the workplace;
- Pay attention to the use of flammable liquids or gases such as disinfectants, anaesthetics, solvents, etc.; the appliance is not designed for use in atmospheres subject to explosion hazards so ensure that any such substances have been completely eliminated or evaporated and dispersed before switching on the appliance (ensure there are no residual substances in the area);
- For the above reason, also pay attention to the presence of endogenous gases and materials impregnated with oxygen (e.g. cotton wadding, gauzes, etc.).

5 - OPERATOR INSTRUCTIONS

5.1 - CONTROLS

The appliance features the following control devices on the front panel (Fig. 6):

- water pump on/off button (1) with orange indicator LED that lights when the pump is running;
- "reverse" button (2) with orange indicator LED; if the selected function requires "normal" rotation direction of tools, this LED must be off;
- the two commands described above are replicated also on the footswitch (Fig. 7), where the green left-hand button (3) controls starting and stopping of the water pump, and the red right-hand button (4) serves to reverse the direction of rotation of the tools;
- "Up ▲" (5) and "Down ▼" keys (6) to edit the field currently selected on the display (8);
- "shift" key (7) to select the field to edit on the display (8);
- two-line display (8); the display shows editable parameters which must first be selected by pressing the "shift" key (the cursor flashes on the selected field) and can then be edited using the "Up ▲" and "Down ▼" keys.
- the footswitch (Fig. 7) is equipped with a yellow button (9) that serves to change the preset program when used in conjunction with the green (3) or red (4) buttons;
it is also equipped with a lever (10) that can be moved to start the motor and control the speed up to the maximum value previously entered in the program.

To change the program with the footswitch first press yellow button (9) and then:

press green button (3) corresponding to Up ▲, if you want to go to the next program;

or press red button (4) corresponding to Down ▼, if you want to go to the previous program.

When you press the yellow button a continuous audible signal sounds to indicate that the green (3) and red (4) buttons are active only in program changeover mode and can therefore not perform their primary function of, respectively, starting/stopping the water pump and reversing the direction of rotation of tools. When you press yellow button (9) again the audible signal stops and the green (3) and red (4) buttons are once more available to perform their primary functions.


Fig. 6


Fig. 7

Figure 8 shows the data that appear on the display as described below:

- 1) sequential program number (from 00 to 31);
- 2) Ncm “xxx” : describes the torque (power) available at the tool tip;
- 3) H2O percentage (%): indicates the quantity of water delivered;
- 4) r1 / r2: indicates which of the two micromotors is selected (r1 or r2);
- 5) transmission ratio between motor and speed reducer or speed multiplier (this parameter is set in accordance with the handpiece installed); reducers are identified by the first digit, multipliers by the second: e.g. 32/1 indicates that speed reducer 32 has been selected, while 1/5 corresponds to speed multiplier 5; when the reducers or multipliers are changed the torque (2) and speed (6) values are updated automatically on the display since these values depend on the currently selected reducer or multiplier;
- 6) RPM “xxxxx” (200 ... 40000 rpm): shows the rpm (speed) delivered by the motor if no reducers or multipliers are selected, or the spindle speed if a reducer or multiplier has been selected.


Fig. 8

5.2 - INSTRUCTIONS FOR USE

5.2.1 - Setting into operation

After performing the installation procedures the appliance can be set into operation:

- 1) equip the micromotor with the required handpiece and tool;
- 2) switch on the appliance by means of the main power switch (6, Fig. 2);
- 3) select the program number by positioning the cursor using the "shift" key (7, Fig. 6) in field 1 shown in figure 8, and then using the "Up ▲" and "Down ▼" keys until obtaining the program required;
- 4) if there is no program saved in the memory that is suitable in relation to the operation you are about to perform, use the edit procedure to change one of the existing programs appropriately (see § 5.2.3);
- 5) push the footswitch lever to start the motor; use the previously described control devices on the footswitch (see § 5.1) to start or stop the water pump and reverse the direction of run of the motor as required.

5.2.2 - Normal deactivation and emergency deactivation procedures

The appliance can be deactivated as follows:

- by releasing the footswitch;
- by setting the main power switch to OFF;
- by unplugging the power cable from the electrical socket outlet.

5.2.3 – Editing a program

To edit a program (see Figs. 6 and 8): position the cursor on the "program n°" field and use the "Up ▲" and "Down ▼" keys to select the number of the first free or unused program.

Once you have selected a program to edit, move around the fields using the "shift" key and change the relative parameter values with the "Up ▲" and "Down ▼" keys; when you have finished the edit procedure to save the program return to the "program n°" field and then quit the program by going to the next or the previous program by pressing "Up ▲" or "Down ▼". If the appliance is switched off before this final step the changes to the parameter values will not be saved in the memory and the new program settings will be lost.

A maximum of 32 programs can be saved in the memory.

ALWAYS SWITCH OFF THE APPLIANCE WHEN YOU FINISH WORKING WITH IT AND FOLLOWING ALL STOPPAGES, IRRESPECTIVE OF THE MOTIVE.

IF YOU NOTICE MALFUNCTIONS, STOP THE APPLIANCE IMMEDIATELY AND DISCONNECT IT FROM THE ELECTRICAL POWER SUPPLY. BEFORE

RECONNECTING THE POWER SUPPLY FIND THE CAUSE OF THE FAULT AND REMEDY IT. IN THE CASE OF PROBLEMS, CONSULT CARLO DE GIORGI.

6 - MAINTENANCE

Periodically check the condition and operation of the accessories and parts subject to wear and tear; in particular check the condition of the burs and the micromotors in order to identify any possible damage.

Periodically check the condition of cable insulation. The appliance cables must be replaced exclusively by suitably skilled and authorised personnel.

6.1 - CLEANING

The appliance does not require any particular maintenance operations, with the exception of regular cleaning of the working areas: before cleaning the appliance always disconnect it from the power supply.

The micromotor supports, stand for saline solution bottle and footswitch IPX7 supplied with the appliance can be cleaned with alcohol;

Any flammable substances utilized for cleaning or disinfecting, or any adhesive solvents, must be allowed to evaporate completely before the appliance can be used.

If you use the tube supplied with the peristaltic pump, flush it with clean water to remove any traces of sodium chloride before stowing the appliance away after use.

6.2 - STERILIZATION

The tube supplied in the peristaltic pump can be cleaned with alcohol or non-toxic sterilizing liquids.

Particular attention must be paid to cleaning and sterilizing the micromotor complete with hose and support.

The micromotor can be sterilized in an autoclave at 135 °C, while the support can be sterilized by autoclaving at 123 °C .

IMPORTANT:

- before placing the micromotor in the autoclave fit the specially designed endcap in the front of the motor (Fig. 9-1 and Fig. 9-2):


Fig. 9-1

Fig. 9-2

- do not sterilize the micromotor with liquids that conduct electricity.

7 - DIAGNOSTICS

If the appliance fails to operate first check the fuses at the rear (7, Fig. 2). If the fuses must be changed refer to the specifications in the table in chapter 2 of this manual. Before checking the fuses, set the main switch to 0 (OFF) and unplug the mains power cable; check the condition of the fuses and replace if necessary. Once the fuse caps have been refitted you can reconnect the mains power cable.

If malfunctions and faulty operation persist, unplug the appliance and consult CARLO DE GIORGI.

WARNING!

**ALL WORK ON THE APPLIANCE MUST BE PERFORMED WITH
THE POWER SUPPLY DISCONNECTED**

8 - ACCESSORIES

8.1 - MOTOR

Autoclavable motor complete with support

8.2 – TUBE FOR IRRIGATION

Alternative tube for the connection between irrigator and handpiece (different with respect to the tube supplied as standard) (See Fig. 10 - A).

8.3 - UNION

Alternative union for the connection between irrigator and handpiece (different with respect to the union supplied as standard) (See Fig. 10 - B).



Fig. 10

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